

TENS Electro  **Medic**

Manage your pain when active or at rest

TRANSCUTANEOUS
ELECTRICAL
NERVE STIMULATOR

**USER'S
MANUAL
IN ENGLISH**

READ INSTRUCTIONS

B E F O R E

OPERATING

Care

1 INTRODUCTION BEFORE USING THE STIMULATOR

1.1 Electro-Medic: A Solid Culture, a Therapeutic Alternative, Devices Conceived from the Finest Technology	3
1.2 Medical context. Usage and Benefits of Neurostimulation	4
1.3 Safety measures. Indications, Contraindications, Precautions, Caution	5-7

2 INTRODUCING THE DEVICE

2.1 Equipment and accessories	8
2.2 Electrical Specifications and Pre-Programmed Settings	9
2.3 Control Pad	10
2.4 Screen Display	11

3 INSTRUCTIONS

3.1 Electrodes. Connection, installation, use, and start-up	12
3.2 Specific Programs. Changing a Program or its Treatment Time, Locking/Unlocking a Program	13
3.3 Adjusting and Locking Intensity Level	14
3.4 Halting or Pausing Stimulation	14
3.5 List of Stimulation Modes	15
3.6 Table of Different Programs	15
3.7 Battery. First Charge, Battery Status, Replacement, Battery Charger	16-19
3.8 Troubleshooting	20

4 DEVICE MAINTENANCE AND CLEANING

21

5 WARRANTY

21

6 DOCUMENT HISTORY

22



BEFORE USING THE STIMULATOR

A CONCEPT ACCESSIBLE TO ALL

INTRODUCTION

1.1

Electro-Medic is the proud Canadian manufacturer of -TENS-, safe and efficient neurostimulator devices. They are compact electrical devices that are therapeutic and portable, offering :



page 3

A SOLID CULTURE IN ELECTROTHERAPY conceived to noticeably reduce acute, chronic, sports-related, post-operative, and post-traumatic pain and much more.

A therapeutic option that is medically recognized and applicable within a clinical environment, at home, and even, depending on the pain to treat, while engaging in everyday activities!

Our models are state-of-the-art, light, compact, and user-friendly and are developed through thorough and consistent collaboration with healthcare providers and their patients to develop a device that is capable of serving an expanded range of physical rehabilitation care.

MANAGE YOUR PAIN WHEN ACTIVE OR AT REST MEDICAL CONTEXT

USING THE TRANSCUTANEOUS ELECTRICAL NERVE STIMULATOR (TENS)

TRANSCUTANEOUS ELECTRICAL NERVE STIMULATION IS SCIENTIFICALLY PROVEN AND
PRESCRIBED BY MEDICAL PROFESSIONALS

Transcutaneous electrical nerve stimulation -TENS- depolarizes peripheral nerve fibres and is transmitted through electrodes placed on the body to reinforce the effectiveness of the natural pain control mechanisms.

TRANSCUTANEOUS ELECTRICAL NERVE STIMULATION APPLIED AT THE SENSORY LEVEL

A tingling sensation is induced to trigger a natural analgesic reaction.

COMPLEMENT OR ALTERNATIVE TO TAKING MEDICATION

Neurostimulation is used for acute postoperative and post-traumatic pain and muscle relaxation and for treating many other types of pain that are not relieved by medication.

ACCESSIBLE TO ALL

The ability to choose the most appropriate TENS program for your type of pain makes TENS one of the most effective non-surgical and non-drug therapy solutions. Transcutaneous electrical nerve stimulation -TENS- is much more than a simple compliment to any traditional treatment. TENS allows you to move more and live better.

1.2
THE
TENS

page
4

1.3 SAFETY MEASURES

INDICATIONS, CONTRAINDICATIONS, PRECAUTIONS, CAUTION

Transcutaneous electrical nerve stimulation meets each user's specific needs. This non-aggressive and drug-free technique is recognized for its lack of side effects under normal use conditions, moderate cost, the small size of the TENS device, and its ease of use. The -TENS- is compact and can be slipped into a pocket, making it easier to carry out daily activities.

BENEFICIAL EFFECTS OF TRANSCUTANEOUS ELECTRICAL NERVE STIMULATION - TENS -

Pain management • Return to an active life • Improved mobility • Accelerated recovery after an accident

TENS INDICATIONS FOR PAIN RELIEF

- Sprain
- Acute or chronic neck, back, and lower back pain
- Tendinitis, epicondylitis, bursitis, capsulitis
- Rheumatism (osteoarthritis, arthritis)
- Tension headache
- Diabetic neuropathy
- Shingles
- Fibromyalgia
- Phantom limb
- Complex regional pain syndrome
- Post stroke or spastic pain
- Multiple sclerosis
- Angina pain
- Childbirth
- Dysmenorrhea
- Intermittent claudication
- Perineal pain
- Interstitial cystitis

- acute, subacute, or chronic pain
- post-traumatic pain
- pre- and postoperative pain

FUNCTIONAL REHABILITATION: better autonomy in performing daily activities

CONTRAINDICATIONS

- Cancer - absolute contraindication: wait until the end of the remission period for metastasis or circulatory cancer (leukemia).
- local contraindication: wait until the end of the remission period for other types of cancer, under interdisciplinary supervision.
- Cardiac pacemaker: absolute or local contraindication, to be used under interdisciplinary supervision with permission from the attending cardiologist.

LOCAL CONTRAINDICATIONS

DO NOT APPLY TO THESE REGIONS

- Transcranial application: the risks of applying the electrodes in a transcranial application are unknown.
- Anterior neck region, carotid sinus (throat): risk of stimulating the vagus or phrenic nerves, pharyngeal muscles, or carotid sinuses.
- Eyes: the risks of treating this part of the body are unknown.
- All types of infections: risk of spreading the infection.
- Skin impairment: psoriasis, eczema, etc.. Caution in the vicinity of open wounds or healing scars.
- Undiagnosed persistent pain.

LOCAL CONTRAINDICATIONS

USE UNDER MEDICAL OR INTERDISCIPLINARY SUPERVISION

- Transthoracic and anterior application of the cardiac region.
- Electronic implant: risk of interference with the implant's function.
- Heart disease: risk that the heart will have difficulty compensating for the high metabolic demand.
- Pregnancy: abdominal and lumbar region during pregnancy. The -TENS- may be used in the lumbar region during childbirth.
- Skin sensory disorder, loss of sensitivity: risk that the patient does not adequately feel the current, which increases the risk of skin burns or irritation.
- Genital organs: requires training.
- DVT/active phlebitis/embolism: risk of thrombus displacement in the bloodstream.
- Hemorrhage: risk of increasing bleeding.
- Epilepsy: local contraindication to the head and neck (cervical region). Precaution for the trunk and limbs. Electrical stimulation could trigger an epileptic seizure.
- Tuberculosis: electrical stimulation in this area may spread the infection during the active phase.
- Pacemaker: absolute or local contraindication, to be used **under interdisciplinary supervision with authorization from the attending cardiologist.**



PRECAUTIONS

- Circulatory dysfunction : stimulation increases the metabolic demand, which may exceed the oxygen supply, thus increasing pain. This may even lead to tissue ischemia or necrosis.
- Skin disease: resistance is decreased, increasing the risk of burns.
- Active epiphyseal plate: risk of impairing bone growth.
- Chest, heart, and lower abdominal area: risk of affecting normal heart function.
- Communication disorder: risk of injury, misunderstanding of use; **under interdisciplinary supervision.**

CAUTION

- Children: keep out of the reach of children.
- Driving a vehicle: never handle the -TENS- device or move the electrodes while driving.
- Risk of skin damage when using carbon electrodes: always use electrode gel.

WARNING

- Do not adjust or perform maintenance on your Electro-Medic -TENS- device while it is in operation.
- Near other equipment. Do not use the device when it is placed near to or above other equipment. If it is necessary to use it in such a configuration, make sure that ALL DEVICES ARE working properly under these conditions.
- Do not use the device at the same time as monitoring equipment (e.g. ECG equipment) operating with electrodes. The signals of the device could interfere with those of the monitoring device.

2

INTRODUCING THE DEVICE

2.1 EQUIPMENT AND ACCESSORIES | MODEL: CARE



THE ELECTRO-MEDIC - TENS- CARE PORTABLE STIMULATOR IS SUPPLIED WITH THE FOLLOWING ACCESSORIES:

- 1 case
- 1 quick start guide
- 1 charger and 1 Li-ion battery
- 2 wires
- 1 pack of reusable square self-adhesive electrodes

CAUTION

Carefully read the instructions for use for the electrodes as explained on the packaging.

We recommend that only accessories authorized by Electro-Medic be used.



2.2

CHARACTERISTICS

ELECTRICAL SPECIFICATIONS

- 2 totally independent stimulation channels
- Constant current for a load up to 1500 Ohms
(Each program has a maximum output amplitude of 38 V)
- Maximum current of one impulse per channel: 0-60 mA
- Impulse duration: 70-250 Us
- Impulse frequency: 1-80 Hz
- Power supply: 1 rechargeable Li-ion battery
3.7 volts/600 mAh
- Weight with battery: 92 g
- Weight without battery: 75 g
- External dimensions: 110 mm (L), 64 mm (W), 11 mm (D)
- Timer: 10-20-30-40-50-60-continuous
- Functional humidity: 5°C to 40°C, 15% and 90% R.H.;
700 hPa and 1060 hPa
- Storage and transport humidity: -10°C to 60°C, 15% and 75% R.H.;
700 hPa and 1060 hPa

TENS SAFETY MEASURES

INTERNALLY POWERED EQUIPMENT, TYPE BF APPLIED PARTS NOT SUITABLE FOR:


- Use with: an anaesthetic mixture flammable with air, oxygen, or nitrous oxide
- Continuous use

CHOICE OF PRESET PROGRAMS

1. Conventional
2. Modulated pulse duration
3. Burst



- 1 INCREASE
Increases the intensity of the left or right channel.
Increase the intensity with caution and as prescribed.
- 2 DECREASE
Decreases the intensity of the left or right channel.
Deactivates lockout.
- 3 ON/OFF BUTTON
Opens and closes the device in addition to halting stimulation at any time.
- 4 PROGRAMS
Choice of 3 preset programs for customized, effective, and safe treatment.
channel 1. See no. 3.5 DIFFERENT POSSIBLE MODULATION MODES
- 5 SELECTION
Confirm and save the timer selection.
- 6 TIMER
Activates the timer to allow users to set the treatment duration.
Choose: 60-minute timer by increments of 10 minutes or in continuous mode **C** depending on medically supervised needs and recommendations.
- 7 PROGRAMS
Choice of 3 preset programs for customized, effective, and safe treatment.
channel 2. See no. 3.5 DIFFERENT POSSIBLE MODULATION MODES
- 8 PAUSE
Pauses the device and returns the intensity to zero.
The device gradually returns to the previous intensity level when the pause button is pressed again.
Timer function will also pause when the device is paused.
- 9 CHANNEL NO. 2 OUTPUT
- 10 CHANNEL NO. 1 OUTPUT




CARE PAD

2.3

page



- 1 PROGRAM LOCKED
Program cannot be changed.
- 2 OPEN CIRCUIT
Electrode disconnected. Other problem likely.
- 3 PROGRAM
Displays the selected program.
 - The left side displays the channel 1 program number
 - The right side displays the channel 2 program number
- 4 INTENSITY
Channel intensity on a scaled bar.
- 5 IMPULSE INTENSITY
Channel intensity represented in numbers.
- 6 LOCKOUT
Indicates if the program is locked.
- 7 TIMER
Displays time remaining.
- 8 DEVICE IN PAUSE MODE
Indicates if the device is in pause mode.
- 9 BATTERY STATUS
Indicates the available current.



SCREEN DISPLAY


Please refer to Nos. 3.2 to 3.8 INSTRUCTIONS, for a detailed explanation of these functions.

3


CONNECTION

3.1 A

CONNECT THE ELECTRODES TO THE WIRES



CONNECT THE WIRES TO THE -TENS- UNIT AT THE INPUT OF EACH CHANNEL



INSTRUCTIONS

For optimal and safe use, use only the accessories provided by Electro-Medic. Furthermore, proper adjustment of the device's intensity and a gradual increase in levels will ensure the comfort, relief, and improvement you desire. We strongly advise you not to increase the levels too quickly.

(Type BF applied part: Electrodes)



B ATTACH ELECTRODES

on non-irritated skin that has been washed with water, cleaned, and dried thoroughly for better electrode adhesion and optimal performance.

Do not position the device where access to the main source of the device would make it difficult to shut the device down quickly if necessary.


D

TURN ON THE DEVICE



To optimize your results, we recommend that you evaluate your pain before, during, and after treatment, using the pain scale

Pain Scale



3.2 SPECIFIC PROGRAMS

CHANGING A PROGRAM OR ITS DURATION

LOCKING/UNLOCKING A PROGRAM

See (3.5) DIFFERENT
POSSIBLE MODULATION
MODES and
(3.6) PROGRAM
SELECTION TABLE

CHANGE A PROGRAM



A competent professional must approve the programmed settings according to your condition.

Press the channel's **PR** button to modify the settings until the desired program is displayed in section no. 3 on the screen or press the up arrow ▲ to start the treatment.

- If you wish to use both channels, select a program in the other channel using the other **PR** button.
- All programs can be used in combination.
- Never change the settings without consulting your therapist.

LOCKING/UNLOCKING A PROGRAM



To lock or unlock a program, press and hold the down arrow on program 2 and the **PR** button of program 2 simultaneously for 10 seconds.

You will see the lock symbol as pictured above in the left section.

- Use the down arrows ▼ to confirm it is working.
- If the lockout symbol **appears**: program locked, cannot be changed.
- If the lockout symbol **does not appear**: program not locked, user can change the program.

CHANGING THE TREATMENT TIME 0-60 MINUTES



Press the Timer button, and the timer will flash. • For the six-step, 10 minute each, treatment option, select 10, 20, 30, 40, 50, or 60 minutes.

The device will automatically count down the time and stop when the time runs out

- For the no interruption, continuous treatment option, keep pressing until you receive the **C** signal.

You will have to stop the device yourself when you consider the treatment time sufficient.

Confirm the chosen option by pressing the **S** button to save or press the up arrow ▲ to start the treatment.

You can repeat the treatment time as needed to experience the desired benefits.

3.3 INTENSITY LEVEL SETTING AND LOCKING

ADJUSTMENT

If you press and hold the up arrow ▲, the intensity level increases faster.

Conversely, if you press and hold the down arrow ▼, the intensity level decreases faster.

Caution: Always increase the intensity level gradually.



LOCKOUT

Depending on the level you want to set, press the up ▲ or down ▼ arrow.

This security lockout symbol indicates that the keypad has been locked out.

The symbol will appear after 10 seconds, once you have selected the intensity level

Therefore, if you wish to increase the intensity, you must first press the down arrow ▼ to deactivate the safety lockout.

You can then increase the intensity level by pressing the up arrow ▲.



3.4 STOPPING OR PAUSING THE STIMULATION

STOP

To stop the stimulation, decrease the intensity with the down arrow ▼ until the intensity returns to 0.00, or press the On/Off button.

PAUSE

You can take a 5-minute break at any time during the treatment.

- If the device is locked, unlock it by pressing the down arrow before pressing Pause
- The timer will stop during the desired pause
- Press Pause again to resume treatment



3.5

DIFFERENT POSSIBLE MODULATION MODES

C conventional STIMULATION

In the conventional stimulation mode, the impulses are continuous with a fixed rate and width.

B BURST

In Burst mode, there are two impulse trains per second (2 Hz).

MW Modulated Pulse duration

The impulse duration decreases from 150 Us to 70 Us in 9 equal stages, then returns to the selected value in 9 stages.

The total cycle time is 3 seconds.

3.6

TABLE OF DIFFERENT PROGRAMS

PROGRAM	PULSE FREQUENCY (Hz)	PULSE DURATION (Us)	TIMER (MIN)	MODE	FIXED
P1 Conventional	80	150	30	C	yes
P2 Modulated pulse duration	80	70-150	30	MW	yes
P3 Burst	2	250	20	B	yes

LEXICON

Hz:	Pulse frequency setting
Us:	Pulse duration setting
Us no. 1:	no.1 Pulse duration setting
Us no. 2:	no.2 Pulse duration setting
Timer:	Timer duration
Mode:	Program modulation
Fixed:	Indicates if the program parameters are fixed or may be changed

BATTERY

3.7

BATTERY STATUS

As displayed:

1/3 of the battery 

2/3 of the battery 

3/3 of the battery 

SERVICE LIFE


The typical service life of a Li-ion battery is approximately:

- Three (3) years or
 - 300 charging cycles¹
- Whichever comes first.

Expect a service life of three (3) years for batteries that do not complete the charging cycle process.


¹ A charge cycle represents a complete discharge followed by a complete recharge of the battery.

BATTERY REPLACEMENT





LI-ION BATTERY CHARGER

LED INDICATOR SPECIFICATIONS

- Adaptor (Model: MDA534627)
- Input: 100 V-240 V~, 50/60 Hz, 150 mA
- Output: 4.2 V , 650 mA
- Red indicator: currently recharging
- Green indicator: no battery or fully charged

CHARGER

- Input: 4.2 V 
- Output: 4.2 V 
- Green indicator: fully charged
- Green indicator: device plugged in without battery

THE PATIENT IS THE INTENDED OPERATOR

The patient can operate the buttons and change the Li-ion battery under normal conditions, as well as maintain the device and its accessories according to the user guide

CAUTION!

- Use only Electro-Medic rechargeable Li-ion batteries
- NEVER reverse the (+) and (-) terminals when connecting, and avoid any contact between the Li-ion batteries and metal objects (necklaces, hairpins, etc.)
- NEVER charge Li-ion batteries for more than 72 hours
- The Li-ion battery charger must comply with IEC 60601-1standards

SAFETY MEASURES

- Do not expose equipment to flame, sunlight, or any other heat source, as this could result in a risk of burns, explosion, or even the release of toxic gases
- Do not store or transport the device with metal objects
- Do not disassemble or modify the device's components
- Avoid contact with water or other liquids

INSTRUCTIONS FOR USING THE CHARGER

- Insert a Li-ion battery. Align the (+) and (-) terminals correctly
- Plug the charger into a standard wall outlet
- A red or yellow LED indicator indicates charging
- When charging is complete, the LED indicator turns green.
Unplug the charger and remove the Li-ion battery.

Indicator colours and their meaning

red	An immediate response from the operator is required
yellow	A quick response from the operator is required
green	Ready to use
Other	Meaning other than these

TIME TO RECHARGE

- A Li-ion battery requires approximately 3.5 hours to recharge.

LI-ION BATTERY

- Voltage limited to 4.2 V
- Rechargeable 3.7 V/600 mAh Li-ion battery

ADAPTOR

- 100-240 V 50-60 Hz, 1.2 A

The adaptor is 2MOPP equipment under IEC 60601-1-1. Equipment approval is valid if used in combination with the adaptor supplied with this equipment.

CAUTION

- This equipment must not be connected with an adaptor other than that supplied with the Electro-Medic equipment.

SAFETY MEASURES

- Do not short-circuit the device
- Do not expose the device to high temperatures
- Use only the specifically recommended charger

SAFETY CLASSIFICATION OF ELECTRO-MEDIC EQUIPMENT

Protection against electrical overload	Electro-Medic device internal power supply
Applied part	Type BF applied part: Electrodes
Protection of electrical equipment against water or fine particle infiltration that may have a harmful effect on the equipment.	IP22
Instructions	Continuous operation

Note: Not intended to be sterilized.
Do not use the device in an oxygen-rich environment.

Do not position the device where access to the main source of the device would prevent the device from being shut down quickly if necessary.

LI-ION BATTERY

- Voltage limited to 4.2 V
- Rechargeable 3.7 V/600 mAh Li-ion battery

SAFETY MEASURES

- Do not short-circuit the device
- Do not expose the device to high temperatures
- Use only the specifically recommended charger

CLASSIFICATIONS

CLASS II DEVICE WITHOUT TYPE BF
APPLIED PARTS NOT SUITABLE FOR:

- usage in the presence of an anaesthetic mixture flammable with air, oxygen, or nitrous oxide
- continuous operation




ADAPTOR

- 100-240 V 50/60 Hz, 1.2 A

The adaptor is 2MOPP equipment under 60601-1-1 Equipment approval is valid if used in combination with the adaptor supplied with this equipment.

CAUTION

- This equipment must not be connected with an adaptor other than that supplied with the Electro-Medic equipment
- Plug the charger into a standard wall outlet, do not use a power bar

PROBLEM	SOLUTION	ACTION
Broken Screen Display problem Replacement part	None	Contact your distributor
The device does not turn on at all	Is the Li-ion battery inserted correctly?	<ul style="list-style-type: none"> • Try changing the Li-ion battery • Charge the Li-ion battery
The sign  appears, intensity increases, but no current	Check the condition of the wires	<ul style="list-style-type: none"> • Try changing the wire • Try the other channel with the same wire
The sign  appears, intensity remains at 1 and does not increase	Check the condition of the wires	<ul style="list-style-type: none"> • Try changing the wire • Try the other channel with the same wire • See if the wire is twisted • See if the connector is damaged
The sign  appears	<ul style="list-style-type: none"> • Are the electrodes on the skin? • Are the electrodes at the end of their service life? • Is contact with the skin good? • Check the condition of the wires • Too much resistance between the electrode and the skin 	<ul style="list-style-type: none"> • Try with carbon electrodes • Change the self-adhesive electrodes • Perform the test on the other channel
Device opens and closes	Check the Li-ion battery or Li-ion battery compartment	<ul style="list-style-type: none"> • Try changing the Li-ion battery • See if the Li-ion battery fits properly in the compartment
Current is unstable	<ul style="list-style-type: none"> • Does the program allow you to feel the current? • Is the wire knotted? • Is the electrode well bonded to the skin? 	<ul style="list-style-type: none"> • Perform the test with P1 • Perform the test on the other channel • Perform the test with carbon electrodes
The stimulation effect is weak or nonexistent	<ul style="list-style-type: none"> • Check the Li-ion battery status • Check the program used • Perform the test on a healthy muscle 	<ul style="list-style-type: none"> • If the problem persists, consult your healthcare professional
The stimulation does not produce the usual sensation	<ul style="list-style-type: none"> • Check the settings • Change the position of your electrodes slightly 	<ul style="list-style-type: none"> • If the problem persists, consult your healthcare professional
The stimulation causes discomfort	<ul style="list-style-type: none"> • The skin is irritated • The contact of the electrode on the skin is not satisfactory • The self-adhesive electrodes are worn • There is not enough electrode gel on the carbon electrodes • The electrode positioning is not optimal 	<ul style="list-style-type: none"> • If the problem persists, consult your healthcare professional

NB: To reduce the risk to the patient and to avoid damage to the stimulator, never increase the intensity (amplitude) to more than 20 mA when testing for possible wire breakage. For more information, do not hesitate to contact your authorized distributor.

MAINTENANCE

DEVICE

Maintenance and cleaning of the device is relatively simple.

- Always store the stimulator and its accessories in the case provided for this purpose
- Never expose the stimulator to water.
- Wipe the device with a damp cloth if necessary.

ELECTRODES

When signs of damage appear, replace the electrodes as recommended in section (3.1) INSTRUCTIONS

a) Self-adhesive electrodes: when the electrodes no longer adhere well to the skin, they can be re-moistened with a few drops of water, then placed back on their plastic film and stored in their airtight bag until the next use.

b) Carbon electrodes: rinse the electrodes when treatment is complete. Never use soap or cleaning products as they may shorten the service life of the electrodes.

WIRES

Disconnect the wires and connections carefully and properly.

BATTERY

To maintain the battery's amperage at its optimum level, we recommend you remove the battery from the device when not used for an extended period of time (approx. 1 month or more).

4

5

WARRANTY

Electro-Medic, the manufacturer, certifies that the -TENS-CARE product is free of material and manufacturing defects at the time of delivery.

Electro-Medic provides a 3 (three) years warranty on all its -TENS- devices, applicable as of the date of purchase of the device.

THE ELECTRO-MEDIC WARRANTY

- applies only to the device.
- does not cover any accessories. (wires, batteries, chargers), which are guaranteed for 3 (three) months.
- can only be claimed by the purchaser of the new product and upon presentation of the proof of purchase.

After verification of the defective device, Electro-Medic will replace the product if it is still under warranty.

The replacement device is covered for the warranty period of the original device.

Any modification, abuse, misuse, or accidental damage or repairs made by a third party will void this warranty.

In the event of a problem and for warranty purposes, the defective device under warranty will be shipped to the sales outlet, a representative of which will follow up with the manufacturer as soon as possible.

6

page 22

DOCUMENT HISTORY

CONCEIVED BY
ELECTRO-MEDIC



650 Industriel blvd suite 100 Blainville Qc Canada J7C 5Y7
1 855 230 6334
info@electromedic.ca

For any other information regarding the use of your Electro-Medic
-TENS- device, please contact your authorized distributor.

INFORMATION ABOUT ELECTROMAGNETIC COMPATIBILITY

Tested and approved according to the safety standards
IEC 60601-1 / IEC 60601-1-2 / IEC 60601-2-10

Wireless communication devices such as cordless home devices,
cellphones, cordless phones and their bases, and walkie-talkies may
affect the operation of the equipment and should be kept at a minimum
distance of 3.3 m

(Note: As shown in Table 6 of IEC 60601-1-2:2007 for ELECTRO-MEDIC
equipment, a typical cellphone with an output of 2 W d=3.3 m with an
immunity level of 3 V/m).

The complete EMC tables are available from tenscare upon request.

The CARE is designed to withstand foreseeable disturbances from
electrostatic discharges (ESD), magnetic fields from the main power
supply, and radio frequency transmitters/receivers such as cellphones.

LEGEND



Read the instruction manual before using the stimulator



Type BF Equipment - Type BF insulated (floating) unit



Dispose of the device, batteries, and accessories in compliance with applicable recycling standards



Class II Equipment



Alternating current



Direct current



IP Rating



Barcode



Serial number



Manufacturer

Manual
version 5.0

Software
Version 5.0

page



